

To copy psychtoolbox across computers:
(tested using Matlab versions 2007a and 2007b, and psychtoolbox version 2.54)

1. Copy psychtoolbox directory onto removable storage. Psychtoolbox folder can be found in the 'toolbox' folder in the MATLAB directory.
2. Paste psychtoolbox folder into toolbox folder on second computer.
3. Once psychtoolbox files are on second computer, in Matlab go to 'File' menu and select 'Set Path'.
4. Click the 'Add folder with subfolders' button, and then direct Matlab to psychtoolbox folder.